
RICHARD WOLF

spirit of excellence



Sterilization Basket System


Sterilization Basket System

The optimum solution

The universal Sterilization Basket System from Richard Wolf provides safe and gentle accommodation for the valuable instrument set during sterilization, transport, and storage.

The universal basket concept allows virtually any basket configuration to be put together.








Top-heavy instruments are securely fixed with an open and a closed silicone holder.

The patented double frame prevents the basket mesh from "becoming frayed" and hence causing injury to the user and damage to the sterilization packaging.

Our double frame therefore ensures maximum safety and a long service life.








Minimal contact surface as a result of the special design of the silicone mounts.


Sterilization Basket System

Easy handling and configuration

System view with sterilization basket, high – Standard (8584.1202) and sterilization basket, high – long (8584.1302)



System view with isterilization basket, low – Standard (8584.1212) and sterilization basket, low – long (8584.1312)








The silicone holders (blue) are based on the instrument diameter and should be ordered separately.

Sterilization Basket System



Ordering information


RICHARD WOLF
spirit of excellence

Picture	Designation	Type
	**Sterilization basket, high - Standard** Dimensions (L x W x H): **530 x 250 x 200 mm** consisting of: Basket lower part, instrument rack high and low, 4 support bars, basket lid high with spring lock and 2 handles	8584.1202
	**Sterilization basket, low – Standard** Dimensions (L x W x H): **530 x 250 x 108 mm** consisting of: Basket lower part, instrument rack low, 2 support bars, basket lid flat with spring lock and 2 handles	8584.1212
	**Sterilization basket, high – long** Dimensions (L x W x H): **725 x 275 x 200 mm** consisting of: Basket lower part, instrument rack low, 4 support bars, basket lid high with spring lock and 2 handles	8584.1302
	**Sterilization basket, low – long** Dimensions (L x W x H): **725 x 275 x 108 mm** consisting of: Basket lower part, instrument rack flat, 2 support bars, basket lid flat with spring lock and 2 handles	8584.1312
	**Sterilization basket, low** Dimensions (L x W x H): **530 x 250 x 108 mm** consisting of: Basket lower part, instrument rack low, 2 support bars, basket lid flat with spring lock and 2 handles, silicone mat "Large Diamond"	85844001

Sterilization Basket System

Accessories

Picture	Designation	Type
	**Support bar** for silicone holder, double	8584.8203
	**Silicone holder, closed, triple** for instruments with ø 8 mm – 14.5 mm	8584.9201
	**Silicone holder, open, triple** for instruments with ø 8 mm – 14.5 mm	8584.9211
	**Silicone holder, closed, for up to 5 items** for instruments with ø 1 mm - 7.5 mm	8584.9501
	**Silicone holder, open, for up to 5 items** for instruments with ø 1 mm - 7.5 mm	8584.9511
	**Silicone holder, closed, for up to 7 items** for instruments with ø 1 mm - 7.5 mm	8584.9701
	**Silicone holder, open, for up to 7 items** for instruments with ø 1 mm - 7.5 mm	8584.9711
	**Silicone stop plate** Ensures the axial position of instruments with smaller diameters	85849005
	Starting Set **silicone holder for urology / arthroscopy** consisting of: 3x 8584.9501, 3x 8584.9511, 1x 8584.9701, 1x 8584.9711	8584.9001
	Starting Set **silicone holder for laparoscopy / surgery / gynecology** consisting of: 3x 8584.9501, 3x 8584.9511, 1x 8584.9201, 1x 8584.9211	8584.9002

Picture	Designation	Type
	**System small part basket** Dimensions (w x h x d): 121 x 35 x 121 mm	8584.3003
	**Locking pin** e.g. for fixation of handles Pack of 2	8584.8205
	**Set of fixation brackets** consisting of: O-ring set, 2x angle bracket, 1x allen wrench	85848210
	**O-ring set** (replacement part) consisting of: 1x O-ring 90 x 3 mm, 1x O-ring 116 x 3 mm	74002311
	**Silicone mat "Large Diamond"** Dimensions (w x h x d): 500 x 18 x 226 mm	85843010

Reprocessing cycle

RICHARD WOLF
spirit of excellence



Specifications subject to change without notice.

J 609.I.18.en.5

Richard Wolf GmbH
Pforzheimer Straße 32
75438 Knittlingen, Germany

Tel.: +49 7043 35-0
Fax: +49 7043 35-4300
info@richard-wolf.com

Co-CEOs:
Jürgen Pfab
Jürgen Steinbeck

Trade Register: Mannheim HRB 510031
VAT no.: DE144521586
Tax no.: 48020/00171

find us on
Social Media
richard-wolf.com